## Calcium

Calcium was performed on a Hitachi Model 737 multichannel analyzer (Boehringer Mannheim Diagnostics, Indianapolis, IN).

Calcium reacts with o-cresolphthalein complexone in the presence of 8 hydroxyquinoline-5-sulfonic acid to form a purple complex. The intensity of the final reaction color is proportional to the amount of calcium in the specimen.

-from Laboratory Procedures Used for the Third National Health and Nutrition Examination Survey (NHANES III) 1988-1994 Elaine W. Gunter, Brenda G. Lewis, and Sharon M. Koncikowski, 1996